

# Intraoperative Imaging at Your Fingertips

# 30%

of pancreatic surgical  
interventions are changed when  
using intraoperative ultrasound.\*

Intraoperative ultrasound guides  
pancreatic interventions by:\*

- Finding tumors not previously visualized by pre-OP imaging.
- Determining presence of true tumor invasions into vessels.
- Visualizing variation in liver blood supply.



Dilated pancreatic duct  
I12C5 transducer



Pancreatic tumor  
I12C5 transducer

## Pancreatic Procedures

- Pancreatic duodenectomy
- Central pancreatectomy
- Distal pancreatectomy
- Pancreatic enucleation

# Unlock Greater Access During Surgical Procedures

## I12C5\*\* MINI-T Transducer

- Flexible imaging for many approaches, including open and hybrid techniques.
- Easy-to-maneuver device with fingertip control.
- Optimal organ contact without excessive pressure.
- Integrated needle guidance for predictable intervention.



\*\* I12C5 is compatible with bk5000 software version 5.16.2 and higher.  
I12C5 has not been licensed by Health Canada.

BK Medical  
8 Centennial Drive  
Peabody, MA 01960  
USA  
T +1 978 326 1300  
bkmedical.com

USA  
Sales & Service  
BK Medical  
8 Centennial Drive  
Peabody, MA 01960  
USA  
T +1 978 326 1300  
F +1 978 326 1399  
bkmedical.com

Europe and Rest of World  
Sales, Service & Design Center  
BK Medical  
Mileparken 34, 2730  
Herlev, Denmark  
T +45 4452 8100  
F +45 4452 8199  
bkmedical.com

**bk**   
medical